
INSTRUCTIONS FOR USE
KYPHON Xpander™ II Inflatable Bone Tamps

INDICATIONS FOR USE
The KYPHON Xpander™ II Inflatable Bone Tamp is intended to be used as a conventional bone tamp for the reduction of fractures and/or creation of a void in cancellous bone in the spine (including use during balloon kyphoplasty with a PMMA-based bone cement that is cleared for use in kyphoplasty procedures) hand, tibia, radius, and calcaneus.

DEVICE DESCRIPTION

Figure 1: KYPHON Xpander™ II Inflatable Bone Tamp

The main components of the Inflatable Bone Tamp (IBT) are a Y-Adapter with the IBT shaft and the IBT balloon, which is located near the distal tip, as shown in Figure 1. The balloon is designed to compress cancellous bone and/or more cortical bone as it inflates. Radiopaque markers located at the distal and proximal end of the deflated balloon allow fluoroscopic visualization of the deflated IBT during positioning.

The shaft contains two non-communicating lumens. The outer lumen is used for IBT inflation, and the central lumen contains a removable stiffening stylet for IBT insertion. An exit marker that is printed on the shaft serves as an indicator of the insertion depth.

The proximal end of the shaft has a Y-Adapter with two ports, the stylet port and the inflation port. The inflation port of the adapter is contiguous with the outer lumen. The stylet port is contiguous with the inner lumen and contains the stiffening stylet. The stylet has a cap at its proximal end for optional stylet removal.

The color coded band (i.e. a yellow band for 10/3 size, a black band for 15/3 size, and a clear/white band for 20/3 size), located just distal to the Y-Adapter, states the balloon size.

CONTRAINDICATIONS

- Bleeding disorder or treatment that increases the chance of excessive bleeding.
- In the presence of active or incompletely treated infection at the site where the bone cement is to be applied.
- Any known severe allergy to contrast media and bone cements.
- Should not be used in vertebral dimensions or fracture pattern do not allow safe placement and inflation of the balloon.
- If the patient has a confirmed pregnancy.
- Instability of posterior wall and/or pedicles.

USA

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

WARNINGS

- Breakage of the device may require intervention or retrieval.
- Do not use this product after the expiration (Use by) date printed on the package. The device may not be safe or effective beyond its expiration (Use by) date.
- For a transpedicular approach, if the pedicle is not large enough or stable enough to withstand the procedure, pedicle fracture may occur.
- Complications that may occur during a parapedicular approach include pneumothorax and bleeding.
- The balloon of the IBT may fail due to contact with bone splinters, bone cement and/or surgical tools.
- Do not inflate the balloon until it has been correctly positioned in the vertebral body. Inflating the balloon prior to full deployment may result in premature balloon failure due to contact between the balloon and the access cannula.
- Inflating the IBT beyond the maximum inflation volume may cause the balloon to rupture before reaching the maximum inflation pressure.
- Inflating the IBT beyond the maximum inflation pressure may cause the balloon to rupture before reaching the maximum inflation volume.
- Never use air or gaseous materials to inflate the balloon.

PRECAUTIONS

- It is important to read the Instructions For Use carefully prior to use.
- Do not use if the package is opened or damaged because product integrity or sterility may be compromised.
- Do not use damaged products. Prior to use, inspect the packaging and product to verify that no damage has occurred.
- Prior to use, the IBT should be examined to verify functionality and ensure that its size is suitable for the specific procedure for which it is to be used.
- The KYPHON Xpander™ II Inflatable Bone Tamps are single use devices intended to contact body tissues. Do not reuse, reprocess, or resterilize. Reusing these devices carries the risk of contamination and may cause patient infection or cross-infection, regardless of the cleaning and resterilization methods. There is also an increased risk of the deterioration of the device performance due to the reprocessing steps, which may lead to patient injury or death.
- Reconditioning, refurbishing, or repair, modification of the product to enable further use is expressly prohibited.
- Do not use this product if you have not been properly trained. The IBT should only be used by physicians who are trained in the techniques of bone tamp use. Physicians using the devices should be familiar with the physiology and pathology of the selected anatomy.
- The IBT should be manipulated only while under fluoroscopic observation with radiographic equipment that provides high quality images.
- KYPHON Xpander™ II Inflatable Bone Tamps must be used with KYPHON® Osteo Introducer® Devices containing the cannula with purple hub.
- The IBT should only be used with a suitable inflation syringe that provides adequate volume control required for inflation of the KYPHON Xpander™ II IBT and is capable of attaining an inflation pressure of 700 psi (48 atm).
- Only inflate the IBT with liquid contrast media: a 60% solution is recommended. Follow manufacturer's instructions for contrast media indications, usage and cautions.
- Risk of balloon rupture increases if IBT comes in contact with bone cement when balloon volume has not been decreased by 1cc from inflated volume.
- The inflation characteristics of the IBT are altered by inflation inside bone.
- Handle and dispose of this product in accordance with accepted medical practice and applicable laws and regulations.

NOTE: Reference the inflation syringe Manufacturer's Instructions For Use.

ADVERSE EVENTS

- Adverse events potentially associated with use of the IBT include:
- Embolicism of fat, thrombus or other materials resulting in symptomatic pulmonary embolicism or other clinical sequelae.
 - Rupture with fragmentation of the balloon of the IBT resulting in retention of a fragment within the vertebral body.
 - Rupture of the IBT causing contrast media exposure, possibly resulting in an allergic reaction or anaphylaxis.
 - Deep or superficial wound infection.
 - Retropulsed vertebral body bone fragments which may cause injury to the spinal cord or nerve roots resulting in radiculopathy, paresis or paralysis.
 - Bleeding or hematoma.
 - Pneumothorax.
 - Pedicle fracture.

DIRECTIONS FOR USE

PREPARATION OF THE IBT

CAUTION: Contrast media may have different viscosity and precipitation levels that may cause slower inflation and deflation times. For this reason, the use of 60% contrast media is recommended.

NOTE: Scales of syringes have been found to use either the milliliter (ml) unit or the cubic centimeter (cc) unit. 1cc = 1ml.

- Twist cap (1) on stylet port (3) to tighten stylet on IBT. Push plunger all the way into locking syringe. Attach locking syringe to inflation port (4) on IBT.
- Pull locking syringe plunger back. Turn plunger to lock it in position on the distal slot in syringe.
- Detach locking syringe from IBT. Set IBT aside.
- Prepare inflation syringe with contrast media according to inflation syringe Manufacturer's Instructions For Use.
- Attach luer connecting port on inflation syringe flexible tubing to inflation port (4) on IBT. The system is now prepared and ready to use.

USE OF IBT

- Select IBT size based on site and treatment goal. Figure 1 and Table 1 define the inflated diameter (D) and inflated length (L) of the IBTs in 37 °C water at inflation volume increments to the maximum inflation volume.
- These dimensions may vary during product use due to local variation in bone structure.

Table 1 - KYPHON Xpander™ II IBT Inflated Dimensions (in 37 °C water)				
Model Number	Inflated Dimensions			
Size 10/3	Volume	Diameter (max of D1, D2)	Length (L)	
	4 cc	2 cc	13.7 mm	16.3 mm
	Max. Inflation Pressure 700 psi (48 atm)	4 cc	17.3 mm	21.6 mm
Model Number	Inflated Dimensions			
Size 15/3	Volume	Diameter (max of D1, D2)	Length (L)	
	5 cc	2 cc	12.0 mm	19.9 mm
	Max. Inflation Pressure 700 psi (48 atm)	4 cc	15.6 mm	25.0 mm
	5 cc		16.8 mm	27.4 mm
Model Number	Inflated Dimensions			
Size 20/3	Volume	Diameter (max of D1, D2)	Length (L)	
	6 cc	2 cc	11.2 mm	23.0 mm
	Max. Inflation Pressure 700 psi (48 atm)	4 cc	14.5 mm	28.1 mm
		6 cc	16.7 mm	32.6 mm

- NOTE:**
- Inflating the IBT beyond the maximum inflation volume may cause the balloon to rupture before reaching the maximum inflation pressure.
 - Inflating the IBT beyond the maximum inflation pressure may cause the balloon to rupture before reaching the maximum inflation volume.

Table 2 - KYPHON Xpander™ II IBT Inflated Dimensions (in 37 °C water and in 5.0 pcf Sawbone)				
Size	Average Volume cc	Average Pressure psi (atm)	Average Volume cc	Average Pressure psi (atm)
10/3	2	135.9 psi (9.3 atm)	2	295.3 psi (20.1 atm)
	4	129.9 psi (8.8 atm)	4	329.3 psi (22.4 atm)
15/3	2	128.3 psi (8.7 atm)	2	272.0 psi (18.5 atm)
	4	121.4 psi (8.3 atm)	4	313.9 psi (21.4 atm)
	5	120.5 psi (8.2 atm)	5	325.6 psi (22.2 atm)
20/3	2	136.0 psi (9.3 atm)	2	283.2 psi (19.3 atm)
	4	130.3 psi (8.9 atm)	4	325.6 psi (22.2 atm)
	6	127.4 psi (8.7 atm)	6	342.8 psi (23.3 atm)

- NOTE:**
- Volume and pressure data were established in in vitro testing. Due to differences in patient anatomy and bone structure, the data is not representative of clinical use.

IBT INSERTION

- An access channel is required for IBT placement.
- Follow the Instructions for Use for the chosen KYPHON® Bone Access Tools to create an access channel into the bone. The KYPHON Xpander™ II IBTs should only be used with KYPHON® Osteo Introducer® Devices containing the cannula with purple hub.

INSERTION SLEEVE WITHDRAWAL

- Without taking the insertion sleeve off the IBT, twist and slightly push the insertion sleeve distally and then pull it proximally over the shaft to expose the distal tip of the balloon. If the insertion sleeve comes off distally while pushing, slide it back proximally over the shaft. Always keep the insertion sleeve on the shaft of the IBT during the procedure.

NOTE: The distal tip of the deflated IBT has reached the distal end of the cannula when the exit marker band (6) on the outer lumen of the IBT enters the proximal end of the cannula.

- Place the deflated IBT into the access channel and position it under image guidance using the radiopaque markers (8). A gentle twisting motion with the forward push can aid insertion.
- While holding the IBT in place, inflate to 44 psi (3 atm) to secure IBT in position. Remove the stylet (2), if desired.

IBT INFLATION

- Inflate the IBT under continuous image guidance. Use the lateral view to monitor distance from the anterior and posterior cortex. Use the AP view to monitor the lateral cortices.
- If continuous imaging is not used, increase the volume in small increments (0.25 – 0.5 cc). Assess tamp position in lateral and AP views before proceeding to further volume increase.
- Stop when treatment goal is achieved: any part of the IBT inflated length contacts cortical bone or maximum inflation volume and/or maximum inflation pressure is attained. (See Table 1)

IBT REMOVAL

Deflate according to inflation syringe Manufacturer's Instructions For Use. Remove the IBT from the bone through the cannula with a gentle twisting motion. When removing the IBT, pull the insertion sleeve back over the balloon.

- If there is resistance, connect the 30 ml locking syringe, pull the syringe plunger back to the "30 ml" mark to create a vacuum, and resume the IBT removal.
- Confirm entry of the balloon of the IBT into the cannula. If the balloon does not move into the cannula, advance the cannula over the IBT to the proximal radiopaque marker. Following cannula advancement, withdraw the IBT through the cannula. If resistance is noted, remove IBT and cannula simultaneously.

CAUTION: Do not withdraw the IBT unless the balloon is fully deflated. Never withdraw the IBT against resistance. Determine the cause of resistance under fluoroscopy and take the necessary remedial actions.

COMPLETION OF BALLOON KYPHOPLASTY PROCEDURE

Following void creation in a pathological fracture of the vertebral body, KYPHON® Bone Cement may be introduced to fill the void using the KYPHON® Bone Filler Device or KYPHON® Cement Delivery System. Refer to the bone cement Instructions for Use, KYPHON® Bone Filler Device or KYPHON® Cement Delivery System Instructions for Use.

CEMENT RESISTANCE TECHNIQUE

If leaving an inflated IBT inside the vertebral body during cement fill on the contra-lateral side, prior to injection of cement, ensure IBT volume is reduced by 1.0cc. It is recommended that the balloon is not placed in contact with the bone cement for more than 5 minutes.

STERILIZATION

Sterilized with irradiation.

HOW SUPPLIED

The KYPHON Xpander™ II IBTs are supplied sterile in a peel-open package. In the event of damage to the sterile packaging, do not use and notify the manufacturer.

MATERIALS NOT SUPPLIED

- Inflation Syringe
- Locking Syringe
- Contrast Media
- Small Bowl

STORAGE

The IBTs should be stored in their original packaging. Proper care should be taken to ensure that the IBTs will not be damaged. Store in a cool (at max. 25 °C / 77°F), dry place.

LIMITATION OF LIABILITY

MEDTRONIC WILL NOT BE RESPONSIBLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, CONSEQUENTIAL, OR EXEMPLARY DAMAGES RESULTING FROM REUSE OF THE KYPHON XPANDER™ II INFLATABLE BONE TAMP.

IN NO EVENT SHALL MEDTRONIC BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, CONSEQUENTIAL, OR EXEMPLARY DAMAGES ARISING OUT OF OR IN CONNECTION WITH THE KYPHON XPANDER™ II INFLATABLE BONE TAMP BASED UPON BREACH OF CONTRACT (INCLUDING BREACH OF WARRANTY).

REQUESTS FOR INFORMATION

For further information, please contact Customer Service, Medtronic, Sofamor Danek USA, Inc., 1800 Pyramid Place, Memphis, Tennessee 38132, Telephone: 800 933 2635 (in U.S.A.), 901 396 3133 (outside U.S.A.), Fax 901 396 0356.

© 2011 Medtronic Sofamor Danek USA, Inc. All rights reserved

Symbols Definition / Definición de los símbolos / Définitions des symboles / Definições dos símbolos

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